

FIXATION EXPERIMENTS for FRESH BLASTODERMS

Grid Box
54 904

October 23, 1975

- A 1. Bouin's Regular Method
- B 2. Bouins + Paraformaldehyde
- C 3. Modified Bouin's + Glutaraldehyde (2%)
- D 4. Modified Bouin's
- E 5. Regaud's Paraformaldehyde + Glutaraldehyde (2%)
- F 6. Bouin's Regular Method + Glutaraldehyde (2%)
- G 7. Glutaraldehyde (12%)
- H 8. Regaud's Paraformaldehyde
- J 9. Osmium in CCl_4 (80%) - 5 min. fixation
- K 10. " - 1 min. "
- L 11. " - 20 min. "
- M 12. Dalton's Osmium (~~weak osmium~~) - 1 hr. fixation
- N 13. Dalton's Osmium + Uranyl Acetate (2:1) (weak osmium)
- O 14. Dichromate (45 min.) then Dalton's Osmium
- P 15. Glutaraldehyde (3%) (30 min.) then Dalton's Osmium
- Q 16. Dalton's Osmium (30 min.) then Glutaraldehyde (3%)
- R 17. Uranyl Acetate + Dichromate (1:1)(30 min.) then add Osmium (2%) (1:1)
- S 18. Dalton's Osmium + Uranyl Acetate (2:1)

